

## Read what one satisfied client has to say

"My name is Don Gulock and I am 81 years old. As a result of an acoustic neuroma, I have a complete loss of hearing in one ear. I have suffered from Single Sided Deafness for 25 years. As I got older, I developed a hearing loss in my other ear. Recently, I was fitted with a Phonak CROS. It has raised my hearing capabilities to a new level. Not only can I hear better, it has none of the limitations I saw with the previous device. The Phonak CROS system is very comfortable. It is suitable to wear in- and outdoors. In addition, it has enhanced capabilities to connect to accessories, such as the ComPilot and TVLink S, which help me interact and communicate better. The Phonak CROS system is wonderful and I, for one, am very grateful it has been developed."

*Don Gulock – Phonak CROS wearer*



---

### Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology we develop innovations that help people hear, understand and experience more of life's rich soundscape.

**Interact freely. Communicate with confidence. Live without limits. Life is on.**

Your hearing care professional:


[www.phonak.com](http://www.phonak.com)

# Phonak CROS




Discover the world's best solution for one-sided hearing loss



## Discover all around hearing delight

Phonak CROS can be a satisfying solution for you, if you have hearing in one ear only. The CROS system wirelessly sends the sound from the side where you cannot hear to your better hearing ear. This gives you access to all around hearing delight.



### When your better ear hears normally

CROS



A small transmitter, placed behind or inside the ear that cannot benefit from a hearing aid, picks up sound and transfers it wirelessly to a receiver on your normal hearing ear. This will

allow you to hear effortlessly even when sounds come from the side where you have no hearing.

### When you have a hearing loss in your better ear

BiCROS



The small transmitter placed behind or inside the ear that cannot benefit from a hearing aid picks up sound and transfers it wirelessly to your ear fitted with a hearing aid. At the same

time this hearing aid amplifies sounds coming directly to this ear. Your hearing and understanding experience will improve and your communication capabilities will be better than ever.

### Phonak CROS H2O

Together with a WaterResistant Phonak hearing aid, this robust little system gives you the freedom to enjoy your favorite activities with great confidence. Be it around water, on the beach, working out, in the rain or if you simply forget to take your system off to take a shower, there is no need to worry.



### Phonak CROS

If you appreciate a small, stylish design with all the benefits that Phonak CROS can provide then this is an attractive solution for you.



### Phonak CROS 13 and CROS 312

For the ultimate in discreet cosmetics a Phonak CROS in-the-ear solution is for you. These tiny custom-made devices come with a choice of battery sizes so that you can choose to have the benefit of a smaller device or a longer lasting battery.

